Seago
SG-921-EK7

-

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Představení

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User
manual

Basic

Seago SG-921-EK7 is a sonic brush designed for children. It has 22,000 oscillations per minute and runs on pencil batteries.

Contents of the package

Toothbrush, three teeth cleaning heads, stand, instructions.

First use

Insert the battery into the brush. Make sure that the brush head is firmly connected to the base. Apply a reasonable amount of paste to the toothbrush. Attach the toothbrush to the teeth and turn it on. The first press triggers the light, the second the cleaning itself. The toothbrush is adjusted to always point out after 30 seconds that you should change the

cleaning area by briefly interrupting the vibrations. Cleaning is accompanied by LED lighting and sound accompaniment. The total cleaning time is set to 2 minutes. Then the toothbrush turns off.

Insert a battery

Release the bottom of the brush by turning it briefly to the left and slide it out. Use a phillips screwdriver to loosen the cover. Insert two AAA batteries into the brush (not included). The correct location of the contacts is shown on the toothbrush. Restore everything to its original state.

Deployment / replacement of the head

Brush heads can only be fitted to the body in the only way that the brush

bristles are parallel to the button. You can also see the correct deployment on the sales packaging of the product. Turn the warhead to the right and take it off. Put on the spare and tighten to the left.

Safety instructions

The toothbrush is intended only for brushing teeth.

If you have recently had surgery in the oral cavity, talk to your dentist about using your toothbrush.

If your gums don't stop bleeding while you're in use, talk to your dentist about using your toothbrush.

If children use the toothbrush, they should at least be checked at the beginning.

Disposal

The product contains batteries and/or recyclable electrical waste. In order to protect the environment, do not dispose of the product with ordinary waste, but hand it over for recycling to a collection point of electrical waste.

Declaration of conformity

The importer,
Satomar, s.r.o.
(karlova 37, 614 00,
Brno) declares that

the Seago SG-921-EK7 toothbrush meets all the essential requirements of eu directives that apply to it. The full text of the Declaration of Conformity can be downloaded

<http://www.satomar.cz> . Copyright ©

2020 Satomar, s.r.o. All rights reserved.

Appearance and specifications may be changed without notice. Print errors reserved. The most up-to-date version of the manual can always be found www.mobilnipodpora.cz

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Copyright © 2020 Satomar, s. R. O.