Seago
SG-977

-

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Představení

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User
manual

Basic

Seago SG-977 is a sonic toothbrush designed for children. It has 24000 vibrations per minute, IPX7 durability, lighting and battery operation.

Contents of the pack

Toothbrush, three heads for brushing teeth, instructions.

First Use

Place battery in brush. Make sure the toothbrush head is firmly attached to the base. Apply a reasonable amount of paste to the brush. Apply the toothbrush to the teeth and turn it on. The first push of the button turns the toothbrush on, the next push controls the backlight and countdown time.

The brush is set to indicate that you should change the cleaning area after 30 seconds, with a short vibration interruption each time. Total cleaning time is set to 2 minutes. The toothbrush is then switched off.

Battery Insert

Release the bottom of the brush (Seago is shown on it) by turning briefly to the left and extending it. Place one AAA battery in the brush (not included). The correct location of the contacts is shown on the brush (the plush part is on the bottom of the brush). Put the bottom of the brush back on.

Deployment / replacement of brush head

Brush heads can only be applied to the body by the only way that the brush bristles are parallel to the button. You can also see the correct deployment on the sales packaging of the product. Gently place the head on the metal outlet of the body. If the warhead is in the right direction, it will be able to push up to the body of the toothbrush. If it was slanted, it wouldn't push against the body.

Safety instructions

The toothbrush is for tooth brushing only.

If you have recently had oral surgery, discuss brushing with your dentist.

If your gums don't stop bleeding when you use it, talk to your dentist about using a toothbrush.

If children use the toothbrush, they should be checked at least initially.

Disposal

The product contains batteries and/or recyclable electrical waste. To protect the environment, do not dispose of the product with conventional waste, but turn it in for recycling at an electrical waste collection point.

Declaration of Conformity

This importer,
Satomar, s. r. o.
(Karlova 37, 614
00, Brno)

declares that the Seago SG-977
brush meets all the essential
requirements of the EU directives
that apply to it.

The complete text of the
Statement of Conformity can be
downloaded at

<http://www.satomar.cz> .

Copyright © 2020 Satomar, s. r. o.
All rights reserved.

Appearances and specifications
can be changed without notice.
Printing errors reserved. The most
up-to-date version of the manual
can be found at
www.mobilnipodpora.cz

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Copyright © 2020 Satomar, s. r. Ö.